

Butter Triglycerides - (T24 - T56)

Column: **MEGA-1 - 0.32mm, 0.10µm, 3m + 1m Built-In Ret. Gap**

Catalog Code: GAP-1-032-010-3-1

Conditions:

Injection: On-Column.

Detector: FID, 350°C.

Oven Temp.: 60°C, 15°C/min, 340°C.

Carrier Gas: Hydrogen, 30kPa.

Peaks:

- | | |
|----------------|---------|
| 1. T24 | 13. T46 |
| 2. Cholesterol | 14. T48 |
| 3. T26 | 15. T50 |
| 4. T28 | 16. T52 |
| 5. T30 | 17. T54 |
| 6. T32 | 18. T56 |
| 7. T34 | |
| 8. T36 | |
| 9. T38 | |
| 10. T40 | |
| 11. T42 | |
| 12. T44 | |

T = Total number of Carbons

