



Dimethylanilines

Column: **MEGA-WAX - 0.32mm, 0.25 μ m, 30m**

Catalog Code: S-WAX-032-025-30

Conditions:

Injection: Split 250°C, 1.0 μ L.

Detector: FID 250°C.

Oven Temp.: 130°C Isothermal.

Carrier Gas: Nitrogen, 1.8 mL/min.

Peaks:

1. 2,6-dimethylaniline
2. 2,4-dimethylaniline
3. 2,5-dimethylaniline
4. 3,5-dimethylaniline
5. 2,3-dimethylaniline
6. 3,4-dimethylaniline

