



## Drugs of Abuse - Basic Drugs

Column: **MEGA-1 - 0.25mm, 0.25 $\mu$ m, 15m**

Catalog Code: S-1-025-025-15

### Conditions:

Injection: Split 280°C, 1.0 $\mu$ L, 1:50 Split Ratio.

Detector: FID 320°C.

Oven Temp.: 120°C, 25°C/min, 310°C.

Carrier Gas: Hydrogen, 30 kPa.

### Peaks:

1. Caffeine
2. Lidocaine
3. Procaine
4. Cocaine
5. Butylantrakinone (I.S.)
6. Heroin
7. Papaverine
8. Etaverine (I.S.)
9. Narcotine

