

Peaks:

- | | | |
|----------------------|------------------------------------|----------------------|
| 1. Ethanoic acid | 7. Hexanoic acid | 13. 2-heptanoic acid |
| 2. Propanoic acid | 8. 2-pentanoic acid | 14. Decanoic acid |
| 3. Isobutanoic acid | 9. 3-hexanoic acid | 15. Undecanoic acid |
| 4. Butanoic acid | 10. Heptanoic acid | 16. Dodecanoic acid |
| 5. Isopentanoic acid | 11. <i>trans</i> -2-heptanoic acid | 17. Tridecanoic acid |
| 6. Pentanoic acid | 12. Octanoic acid | |

Free Acids

Column: **MEGA-ACID (FFAP) - 0.32mm, 0.25µm, 25m**

Catalog Code: C-ACID-032-025-25

Conditions:

Injection: Split 250°C, 1.0 µL, 1:50 Split Ratio.

Detector: FID 250°C.

Oven Temp.: 80°C, 2.5°C/min, 210°C.

Carrier Gas: Hydrogen, 40 kPa.

